


ZNCL-20L

SHORT PATH DISTILLATION


2L-20L for your choice,
Can meet customer's customized needs



ZNCL-20L

PRODUCT COMPONENTS

- 1. Thermometer Probe
- 2. Distillation Head
- 3. Receiving Bottles
- 4. Cold Trap
- 5. Digital Thermometer
- 6. Boiling Flask
- 7. Heating Mantle
- 8. Iron Stand
- 9. Lifting platform
- 10. Receiving Bottles




MULTIPLE FUNCTION
GOOD PRACTICAL

PRODUCT FEATURES




Borosilicate Glass



Rapid Separation



Reliable Sealing



Fast Heating











Less Residue




Support Customization

MULTIPLE FUNCTION
LED DISPLAY

PRODUCT FEATURES

- | | |
|--|--|
|  Heating
Button |  Stirring
Button |
|  Set
Temperature/sSpeed |  Decrease
the Value |
|  Reduce
the Temperature |  Increase
the Value |
|  Timer
button |  LED
Display Screen |




ZNCL-20L

PRODUCT DETAILS

01 Heating Mantle

Maximum temperature: 350 °C
Magnetic stirring



02 Boiling Flask

Borosilicate glass material
20L capacity



03 Distillation Head

It is connected with cooling chiller,
which can make the material
cool and condense quickly.



04 Thermometer Probe

Thermometer microprocessor
-50~1300°C -58~2372°F



ZNCL-20L

GROUPED EQUIPMENTS

Vacuum Pump



Circulating Water Vacuum Pump



Rotary Vacuum Pump



Circulating Water Vacuum Pump

Temp. Control Equipment



Cooling Chiller






DFY Heater-Chiller



OPTIONAL ITEM
CUSTOMIZED SERVICE

THERMAL JACKET

-  Reduce heat loss
-  Low cost
-  Better effect



MULTIPLE FUNCTION
LED DISPLAY

ALL MODELS DISPLAY



SPD-2L

SPD-5L

SPD-5L



SPD-10L

SPD-20L

SPD-20L