

High Temperature Circulating Water/Oil Bath

GY-5L~100L

5L/10L/20L/30L/50L/80L/100L FOR YOUR CHOICE
NO-LOAD MAX. TEMPERATURE: 200°C

PARAMETERS

PRODUCT INFORMATION



Model	GY-5L	GY-10L	GY-20L	GY-30L	GY-50L	GY-80L	GY-100L
Liquid Storage Volume (L)	5	10	20	30	50	80	100
No-load max. Temp.(°C)	200						
Heating power (kW)	2	2.5	3.5	4	5	8.4/9	8.4/9
Power(V/HZ)	220 50/60			380 50/60 (50L/80L)			
Pressure (MPa)	1.45						
Flow Rate(L/min)	15~20		20~25				
Temp. control range(°C)	Water bath : RT~99 Oil bath : RT~250						
Optional configuration	Explosion-proof motors and electrical						
Liner size(cm)	25*22	25*25	30*30	35*30	45*35	45*45	50*50
Overall dimensions(cm)	60*40*68	60*45*70	63*48*70	70*54*85	75*58*85	86*68*90	90*70*95

ADVANTAGES

SOURCE FACTORY



15-25L/min Flow



2-9kW Heating power



RT~250°C for oil bath



RT~99°C for water bath



Sealed magnetic pump cycles



Closed structure



Smooth operation




Anti-Corrosion

PRODUCT DISPLAY

Optional configuration: Explosion-proof motors and electrical



Front display



Internal structure

PRODUCT DETAILS

Source factory, professional R & D team



01

Water/oil bath



02

Universal Caster



03

Liquid outlet

COMPONENTS

PROFESSIONAL R & D TEAM



Digital Control



Technical Support



Support for Custom



APPLICATION

PROVIDE HEAT SOURCE



GY-5L + Jacketed Glass Reactor



GY-10L + Stainless Steel Reactor



GY-20L + Crystallization reactor



GY-30L + Decarboxylation reactor

CONTROL PANEL

INTELLIGENT CONTROL



Set Button



Run/Stop Button



Reduce Temperature



Increase Temperature



Heat On/Off



Cycle On/Off