

ZZKD

Low Temperature Coolant Circulation Pump

DLSB 5L

5L~100L / -10°C~-120°C FOR YOUR CHOICE

REGULAR AND EXPLOSION-PROOF MODELS ARE AVAILABLE



PARAMETERS

PRODUCT INFORMATION



Model	DLSB 5L-10°C	DLSB 5L-20°C	DLSB 5L-30°C	DLSB 5L-40°C
Storage Volume	5L			
Low temperature range	-10°C~RT	-20°C~RT	-30°C~RT	-40°C~RT
Temperature stability	±2°C			
Refrigeration capacity	280~620W	172~2228W	320~1150W	10~1300W
Rated Power	420W	635W		954W
Rated Voltage/Frequency	220V 50HZ			
Max. Input Current	2.1A	2.9A	4.3A	
Flow	20 L/min			
Lift	4~6 m			

ADVANTAGES

SOURCE FACTORY



Brand Compressor



Fast Cooling



Anti-Corrosion



P.I.D. Control



304 Stainless Steel



Easy to Move



Smooth Operation



Low Noise

CONTROL PANEL

Liquid crystal display, intelligent Control



Set button



Raise the Temperature



Reduce the Temperature



Circulation button



Refrigeration button



Power switch

COMPONENTS

PROFESSIONAL R & D TEAM



Digital Control



Technical Support



Support Customization

Liquid outlet

Liquid inlet

Control panel


Cooling window

Anti-corrosion shell

Universal Caster


PRODUCT DETAILS

Source factory, professional R & D team



01

Cooling coil, high cooling efficiency



02

Radiator window, fast heat dissipation



03

Liquid inlet & Liquid outlet

USAGE SCENARIO

WIDELY USED IN VARIOUS FIELDS



Medicine Industry



Science Experiment



Chemical Industry



Research Institute

APPLICATION

PROVIDE COLD SOURCE

Rotary Evaporator

Jacketed Glass Reactor

Closed Loop Extractor

Centrifuge

