

HEATING, COOLING, CIRCULATION

GDX 20L-30°C

GDX HEATER-CHILLER

✓ 5L~100L

✓ -120°C~200°C



Stainless steel liner

Smart control panel

Precise temperature control

Brand compressor

Closed structure

Safe and energy saving

GDX HEATER-CHILLER

PRODUCT DETAILS

1

Control Panel

Intelligent controller to ensure safer operation.



2

Oil Filler


Stainless steel material, acid and corrosion resistance.



3

Cooling Window

Fast heat dissipation, built-in brand compressor.



GDX HEATER-CHILLER

PRODUCT STRUCTURE



1 Control panel

2 Liquid level line

3 Corrosion-proof shell

4 Universal wheel



5 Oil filler

6 Liquid outlet

7 Liquid inlet

8 Cooling window

GDX HEATER-CHILLER

PRODUCT DISPLAY

5-100L

-120°C~200°C

304

① Explosion-proof GDX

② Conventional GDX



GDX HEATER-CHILLER

PRODUCT PARAMETER



PARAMETER

Model: GDX 20L-30°C
 Temp. Range: -30~199°C
 Cooling capacity: 550~3460W
 Heating power: 3000W
 Total Power: 4300W
 Flow: 25L/min Lift: 6m
 Voltage: 220 V
 Current: 22A

ⓘ Different temperature and capacity can be customized.

Model	GDX-5/10	GDX-5/30	GDX-10/30	GDX-20/30	GDX-30/30	GDX-50/80	GDX-100/60
Temp. Range (°C)	-9~199	-30~199	-30~199	-30~199	-30~199	-80~199	-60~199
Refrigerant	R404	R404	R404	R404	R404	R404 R23	R404 R23
Cooling Capacity (W)	380~1200	150~2425	275~2150	550~3460	850~9200	420~14500	2400~18500
Heating Power (W)	1000	1000	2000	3000	4000	6000	9000
Total Power (W)	1600	1800	2900	4300	6400	16400	22000
Circulating Pump Power (W)	100	100	100	100	100	280	280
Flow (L/min)	20	25	20	25	25	35	35
Lift (m)	6	6	6	6	7	11	11
Voltage (V)	220	220	220	220	220	380	380
Current (A)	9	10	16	22	31	28	46

PROVIDE BOTH COLD AND HEAT SOURCES

APPLICATION



GDX HEATER-CHILLER

BRAND COMPRESSOR



COPPER Material



EFFICIENT Compressor



FASTER Cooling

